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
The software tool OWISoft allows the configuration and operation of several OWIS® control units of the type PS 90 (V 2.4 and higher), PS 35 (V 5.0 and higher), PS 30 (V 2.4 and higher), PS 10 (V 1.7 and higher), DC 500 (V 2.5 and higher), SMS 60 (V 2.4 and higher). It comfortably offers the readout and correlation of the parameters from the detecting chip OWISid together with the PS 10. This further simplifies the configuration of the OWIS® positioning system. The settings can be stored. All linear and rotative OWIS® positioning units are stored in the program database. They can be called, modified and then operated by the control unit. The program can execute and asynchronously approach up to 12 axes. The joystick modus is supported. Furthermore, the PS 90, PS 35, PS 30 and the PS 10 can operate in- and outputs.

Technical Requirements

OWISoft can be installed on a commercially available PC.

Hardware Configuration


At least a Pentium 233 MHz (or compatible) with 256 MB RAM and 200 MB fixed-disk storage is to be used.

Interfaces

The PS 90 and PS 35 are connected via USB or RS232 to the computer. The PS 30 is a plug-in card (PCI bus). The PS 10 is connected via USB to the computer. The USB and PCI driver are available. With the installation of the USB(PCI) driver a virtual serial port is installed. The DC 500 and the SMS 60 can be connected over a RS232 or a GPIB interface (only with the National Instruments drivers – NI-GPIB).

Software

OWISoft is written for Windows operating systems and runs under Win98, WinME, Win2000, WinXP, Vista (32/64 bit) and Win7 (32/64 bit). The SDKs for C, C++ and LabView (V 5.0 and higher) are available.

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The recommended screen resolution is 1024x768 pixels, with 24-bit colour depth (or higher).
The connected controls will be configured and then connected to OWISoft. The configurations can be saved.

Afterwards, the actuators connected to the control units are defined. At the end, the parameters of the individual axes can be set and saved. The parameter sets of OWIS® actuators are stored in a database, they can be loaded and modified. Additionally, it is possible to parametrize unknown positioning units (other manufacturers).

After that, the actuators are ready for use.

Using the two applications Free Positioning and Auxiliary the user can directly call functions or enter adjustment movements. They can not be stored. Furthermore, three storable application interfaces are available. With the application Chain a simple, process oriented program can be set in a table (like EXCEL). Meander is an application used to implement two-dimensional patterns or three-dimensional grids. With the application Program the user can generate, test and process complex motion sequences in a script language.

Additionally, OWISoft offers the possibility to define 11 function keys in order to have quick access to OWISoft commands or external applications.