


Light. Measurement. Excellence.



ULS300

Variable Radiance
Uniform Light Source



www.bentham.co.uk

Achieving high quality images with sensors and cameras requires accurate uniformity analysis at the pixel level, particularly where wide-angle lenses are used. Let the unparalleled spatial uniformity of the ULS300 underpin your imaging system characterisation.

- ▶ Pixel normalisation
- ▶ Flat fielding
- ▶ Distortion correction
- ▶ Defect pixel analysis



Superlative Uniformity

A unique design based on external source and multiple fibre bundle illumination of the front hemisphere ensure ultimate spatial uniformity.

NMI Traceable Calibration

Supplied with spectral radiance calibration, impart NMI traceability to your measurements.

Variable Radiance


Tune the ULS300 to the levels you want with ease using a micrometer-controlled slit.

Constant Spectrum


Using a single source with micrometer slit attenuation, the spectral distribution of the ULS300 does not vary with radiance level.

The design of the ULS300 includes an external light source, bi-lateral slit attenuation and eight branch fibre bundle illumination of the front hemisphere. As radiance is varied, uniformity and spectral distribution remain constant.


1. Superlative Spatial Uniformity



2. Constant Spectral Distribution



Achieving ultra-high spatial uniformity whilst maintaining high radiance levels requires a detailed optical design. With the ULS300 benefit from unparalleled spatial uniformity, continuously variable radiance and constant spectral distribution.



1

Single Source

An ellipsoidal reflector ensures maximised coupling from a stable 150W halogen lamp.

2

Fibre Bundle Delivery

Eight-branch fibre bundle to the front hemisphere ensures uniform light distribution in the integrating sphere and superlative uniformity.

3

Integrating Sphere

Coated with high, diffusely reflective barium sulphate, a 300mm diameter sphere offers uniform illumination, even over a 100mm diameter port.

4

Port

The 100mm port may be reduced to 50mm according to application, resulting in higher radiance levels.

5

Monitor Port

A monitor port DH400_VL photometric detector with ORM400 for real-time reporting of luminance (cd.m-2).

6

SMA Port

An SMA port allows coupling to an array spectrometer for optional spectral monitoring of the source.

7


Radiance Control

Micrometer-actuated, bi-lateral slits located between source and fibre allow for precise radiance control.

8

Radiance Stability

Excellent stability is guaranteed by constant current operation from the PSU_610 stabilised lamp power supply.



Electro-Optical	
Lamp type	Grit-blasted quartz tungsten halogen lamp, G6.35 base
Nominal lamp power and voltage	150W, 24V
Operating current	6.3A
Expected lifetime	2000 hours
Coupling to fibre	Ellipsoidal reflector
Coupling to sphere	8-branch glass fibre
Sphere coating	BaSO ₄
Sphere diameter	300mm
Window diameter	100mm (50mm reducer)
Variable radiance actuator	Micrometer controlled bi-lateral slit
Slit width range	0-10mm

Optical Performance	
Spectral range	350-2000nm
Peak spectral radiance (typ.) 100mm	170 mW·m ⁻² ·sr ⁻¹ ·nm ⁻¹ at 750nm
Peak spectral radiance (typ.) 50mm	285 mW·m ⁻² ·sr ⁻¹ ·nm ⁻¹ at 750nm
Luminance range (typ.) 100mm	0- 8500 cd·m ⁻²
Luminance range (typ.) 50mm	0- 14500 cd·m ⁻²
Correlated colour temperature (typ.)	3300 ±20K over full range
Chromaticity coordinates, CIE 1931 & 1976 (typ.)	x = 0.420, y = 0.410 u' = 0.240, v' = 0.520
Uniformity	0.05% over 100mm diameter aperture, independent of luminance setting

Calibration	
Measurement type	Spectral radiance
Wavelength range	350nm-800nm
Wavelength interval	5nm
Calibration frequency	Recommended after 200 hours burn-time
Traceability	Physicalish Technische Bundesanstalt (PTB, Germany)

ULS300	
Connector	4mm socket
Dimensions	290L x 125W x 110H (mm)
Weight	10 kg

Power Supply	
Output voltage & current	0-10.4A in remote use Max Voltage: 26V
Resolution	0.1A local; 0.001A over SCPI protocol
Max output power	250W
Operating conditions	Ambient temperature: 0°C - +40°C Storage temperature: -20°C - +85 °C
Drift	0.05% of rated output (Over an 8-hour interval with constant line, load and temperature, after a 30-minute warm-up period)
Temperature stability	30ppm/°C
Accuracy	±(0.04% +0.05A)
Output power ripple	<0.1%
Control	Front panel / USB (implementing SCPI protocol)
Front panel control features	Power on/off lamp, select pre-set current, define current
Connector	4mm socket
Display	LCD display of set current, actual current, lamp voltage, burn time and power
Dimensions	450L x 300W x 130H (mm)
Weight	5 kg
Power supply	Mains input 110/220V 50/60Hz

Part Number	Description
ULS300	Variable radiance uniform light source
PSU_610	Current stabilised lamp power supply
ORM400	Dual-channel picoammeter with display
DH400_VL	Precision photometric detector

UK Contact Information

Bentham Instruments Limited
2 Boulton Road
Reading
Berkshire
RG2 0NH
United Kingdom

+44 118 975 1355
sales@bentham.co.uk
www.bentham.co.uk



Distributor Contact Information

AUSTRALIA

Photometric Solutions Intl Pty Ltd
sales@photometricsolutions.com

AUSTRIA

LOT-Quantum Design GmbH
info@lot-qd.de

BELARUS

UE PROFCON
marketing@profcon.by

CHINA

IVT Solar Pte Ltd
sales@ivt-solar.com
Sensing Instruments Co., Ltd,
sales@sensingm.com

FINLAND

Oy Mitaten Finland Ab
sales@mitaten.fi

FRANCE

Trioptics GmbH
tristan.athanaze@trioptics.fr

GERMANY

LOT-Quantum Design GmbH
info@lot-qd.de

INDIA

ATOS Instruments Marketing Services
atos@atosindia.com

INDONESIA

Industrial Vision Technology (S) Pte Ltd
wangyc@visiontec.com.sg

ISRAEL

IL Photonics
moshe@ilphotonics.com

ITALY

2M strumenti
info@2mstrumenti.com

JAPAN

Soma Optics Ltd
sales@somaopt.co.jp

KOREA

Won Woo Systems Co. Ltd
sales@wonwoosystem.co.kr

MALAYSIA

Industrial Vision Technology (S) Pte Ltd
wangyc@visiontec.com.sg

NEW ZEALAND

Photometric Solutions Intl Pty Ltd
sales@photometricsolutions.com

PORTUGAL

Lasing SA
info@lasing.com

RUSSIA

Ultratherm LLC
info@ultratherm.ru

SINGAPORE

Industrial Vision Technology (S) Pte Ltd
wangyc@visiontec.com.sg

SPAIN

Lasing SA
info@lasing.com

SWITZERLAND

LOT-Quantum Design AG
dumouchel@lot-qd.ch

TAIWAN

Asia Qtech Instrument Inc
service@qtechinstrument.com

U.S.A.

Market Tech, Inc.
info@markettechinc.net