

ECGt5050/AI/HL/RT

1008841

Technical Specifications

Size and Dimensions	
footprint; height	50 mm x 50 mm; 17 mm
distance center of rotation to bottom	77 mm
maximum installation space	50 mm x 67.6 mm; 19.4 mm
weight (aluminium version)	137 g
height	17 mm
cable	50 cm cable with connector
Materials	
positioner body	aluminum
actuator	PZT ceramics
connecting wires	copper, jacket: RT: PTFE, HV/UHV: fiberglass
Coarse Positioning Mode	
travel range (step mode)	10 °
maximum drive velocity @ 300 K	3 mm/s
input voltage range	0 - 45 V
Fine Positioning Mode	
fine angular positioning range @ 300 K	1.4 m °
fine positioning resolution	μ°
input DC voltage range	0 - 45 V
Load (@ ambient conditions)	
maximum load	1 kg
maximum dynamic torque around axis	35 Ncm
General Specifications	
environment	/RT
operating temperature range	17 - 25 °C
operating humidity range	40 - 60%

