

ECR5050hs/AL/NUM/RT

1012480

Technical Specifications

Size and Dimensions	
footprint; height	50 mm x 50 mm; 15 mm
weight (aluminium version)	102 g
height	15 mm
cable	50 cm cable with connector
Materials	
positioner body	aluminum
actuator	PZT ceramics
connecting wires	copper, jacket: RT: PTFE, HV/UHV: fiberglass
Coarse Positioning Mode	
travel range (step mode)	360 ° endless
maximum drive velocity @ 300 K	15 °/s
input voltage range	0 - 45 V
Fine Positioning Mode	
fine angular positioning range @ 300 K	15 m°
fine positioning resolution	μ°
input DC voltage range	0 - 45 V
Position Encoder	
readout mechanism	optoelectronic sensor
encoded travel range	360°
sensor resolution	0.01 m°
sensor power (when measuring)	300 mW
repeatability	1 m° (bidirectional)
Load (@ ambient conditions)	
maximum load	2 kg
maximum dynamic torque around axis	1 Ncm
Error Motion	
wobble	< 1 mrad
General Specifications	
environment	/RT
operating temperature range	17 - 25 °C
operating humidity range	40 - 60%

