Dynamic Infrared Twyman-Green Interferometer

Instantaneous Acquisition

The PhaseCam® 6100 SWIR is a highly compact, lightweight dynamic laser interferometer with fully motorized controls for measurement of optics and optical systems. Operating at a wavelength of 1550 nm the PhaseCam 6100 SWIR is ideal for long optical path-length measurement of large focal optical systems such as concave telescope mirrors and lens systems, and for testing small aperture afocal components such as flat mirrors and collimators.

The PhaseCam 6100 incorporates Dynamic Interferometry® technology, using a single camera, high-speed optical phase sensor to make wavefront measurements thousands of times faster than a temporal phase shifting interferometer. Because acquisition time is so short, the PhaseCam can measure under almost any conditions, without vibration isolation, making the PhaseCam ideally suited for the production floor, in clean rooms and in environmental test chambers.

Compact and lightweight, the PhaseCam 6100 was designed with performance and remote measurement in mind. Moving the system to reconfigure a test set is simple and easy, and isolation equipment is not required. Fully motorized controls make it easy to operate the system in remote locations.

The PhaseCam 6100 is a turnkey instrument that includes the interferometer, 4Sight™ advanced wavefront analysis software, and complete computer system. Samples with reflectivity from 1% to 100% can be measured with a simple adjustment. The diffraction-limited, custom-designed optical system maximizes sampling of the full aperture of the test part. The easy to use, vibration insensitive PhaseCam 6100 ensures rapid and accurate data acquisition.

Industry Leading Analysis, Standard

4Sight wavefront analysis software features a user-friendly interface with unmatched simplicity, analysis features and graphical displays. The Measurement Screen aids alignment and execution of single, averaged, burst or continuous data acquisition. The Measurement Flow interface lets you visualize the entire measurement process, from raw acquisition through masking, reference subtraction, terms removal, etc. Zernike, Seidel, geometric and diffraction analyses are easy to perform. Comprehensive data sharing capabilities let you read, write, save and print from most file types, including MetroPro IDL®, MatLab®, Opticode®, Vision®, HDF5® and CodeV®. Generating phase movies to characterize deforming surfaces and moving parts is simple and straightforward.

Complete Measurement System



PhaseCam 6100 SWIR

FEATURES

- Vibration Insensitive Dynamic Operation
- 30 µsec Data Acquisition Time
- 1550 nm Wavelength
- Easy Sample Reflectivity Adjustment
- Outstanding Data Analysis and Visualization Software

APPLICATIONS

- Meter-Class Telescope Optics
- Quality Verification of Optical Components
- Vacuum and Environmental Chamber Testing
- Production Floor Quality Control



Specifications

Configuration Model 6100 SWIR

Description Vibration insensitive dynamic Twyman-Green interferometer

Acquisition Mode Single camera, high-speed optical phase sensor Laser Source 1550 nm

 Laser Source
 1350 mm

 Maximum Output
 < 30mW</td>

 Maximum Cavity Length
 > 60 m

 Beam Diameter
 9 mm FWHM

 Polarization
 Circular

 Field of View
 9 mm

Focus Range ±12.5 mm, optical magnification dependent

Pupil Magnification 1X Fixed

Fringe Contrast
User adjustable for reflectivity from 1–100
Camera
512 x 512 pixels, 10-bit standard
Data Array
User selectable full, half, quarter data arrays

Motorized Controls Focus, reference beam block, aperture block, contrast adjustment

Computer System High performance PC

Operating System Windows® 7

System Software 4Sight™ Analysis Software

Instantaneous Phase Shifting data acquisition

Reference generation, subtraction, data averaging, masking

2D and 3D surface maps

Zernike / Seidel / Slope / Geometric / Fourier Analysis

Fiducial aided data set mapping

HDF4 / HDF5 data format standard, others supported Absolute sphere, prism & corner cube analysis

Multiple sub-aperture analysis

Upgrades free during warranty period

Physical Envelope < 40.1 x 15.7 x 9.1 cm (15.8 x 6.2 x 3.6 in)

 $\begin{array}{ll} \mbox{Weight} & < 4.5 \mbox{ kg (10.0 lbs)} \\ \mbox{Power consumption} & < 750 \mbox{ Watts} \end{array}$

Temperature Range Operational: 16–27° C (60–80° F), non-condensing

Storage: -1-38° C (30-100° F), non-condensing

Options

Divergers Range of lenses from f/1 to f/32

System Performance

Acquisition Rate \geq 15 frames/sec live video

≥ 30 frames/sec acquisition with post processing

Minimum Exposure 30 μsec
Sample Reflectivity 1–100%
RMS Repeatability < 0.001 wave**
RMS Precision < 0.002 wave***

Warranty One Year, limited, on-site system installation and operator training



** One sigma for RMS of 10 data sets of calibration mirror, each data set being an average of 16 measurements.

****Average RMS of the pixel by pixel difference of 10 data sets between measured surface and the calibrated surface. Each data set is an average of 16 measurements. Calibrated surface is the average of all 160 measurements.

Patent 7,230,717. Other patents may apply.

 $\label{phaseCam} \mbox{ and Dynamic Interferometry are registered trademarks of 4D Technology Corporation.}$

MetroPro IDL, MatLab,Opticode, Vision, HDF5, CodeV and Windows are registered trademarks of their respective owners.

All specifications subject to change without notice.

1.32

→ 1.01 -



14.75

21

6.15

 $2.56 \rightarrow$

(

←1.98>

INVISIBLE LASER RADIATION, AVOID EXPOSURE TO BEAM. CLASS 3B LASER PRODUCT <30mW at 1550nm

11.14.2017 ■ ©4D Technology Corporation