

## A Unique Time-of-Flight Micro-Optical Velocimeter

The microV is about the size of a pen and requires no calibration or alignment. It was designed for applications where close-range flow measurement and precision positioning are required. The sensor can operate in environments where hot wire anemometers are difficult or impossible to use, such as two-phase flows, temperature varying flows, low speed flows, and high vibration.

## ADVANTAGES OF THE MICROV:

- Reliable, accurate, rugged
- Extremely compact
- No alignment needed
- No calibration needed
- Makes accurate measurement of fluids of varying temperature, pressure, and density
- Can measure flow speed and flow direction
- Computer controlled 1-D traversing system
- Automated profile measurement
- Battery operated option
- Waterproof and temperature resistant housing

## **APPLICATIONS INCLUDE:**

- Sensors embedded into flow models for precise positioning
- Fluid mechanics, turbulence, oceanography, and atmosphere studies
- Micro channels
- Wind, water, and oil tunnels and channels
- Speed of boats, ships, sea bearing vessels



The microV is the smallest non-intrusive flow velocity sensor available in the industry.



The microV System consists of a microV probe, driver electronics, and BP-microV acquisition hardware and processing software (computer and traverses are optional).



The microV probe is easy to incorporate into flow models with its small size, simple geometry, and rugged construction.

MEASUREMENT SPECIFICATIONS		
Velocity range	-100 to 100 m/sec	
Repeatability	99.8%	
Accuracy	99.7% typical	

PROBE VOLUME	
PV dimensions (x by y by z)	200 x 100 x 300 μm
Standoff distance (air/water)	7.4 / 10 mm 0.29 / 0.39 inches

PROBE SPECIFICATIONS	
Probe weight	25 g
Dimensions	9 (dia) x 25.4 mm 0.35 (dia) x 1 inch

LASER SPECIFICATI NS	
Laser power	110 mW
Wavelength	658 nm
Laser type	Class IIIb

OPERATING PARAMETERS		
Temperature	0 to 55°C	
Pressure	Up to 1 bar	
PC requirements	Laptop or PC	

## **OPTIONAL FEATURES**

Micro-traversing stage for profile measure-

Integrated system housing (for flow model integration)

Steel-jacketed cable

POWER SUPPLY 12 VDC Universal

One or more of the following U.S. Patents apply: No. 6,654,102, 6,580,503, 6,608,6668, 6,717,172, 6,956,230 V201511





SCIENCE